



Weitere Produkte  
auf Anfrage!



### Transurethraler Katheter Tiemannspitze

Tiemannspitze, Länge 40 cm  
Ballonvolumen 10 – 30 ml

- großes Lumen
- weicher Ballon
- reines Silikon
- bunte dilatierende Spitze
- große Augen
- Röntgenkontraststreifen

REF	Größe	Ballon	PZN
01/01/12/12/10	CH 12	10 ml	01150606
01/01/12/14/30	CH 14	30 ml	01150629
01/01/12/16/30	CH 16	30 ml	01150635
01/01/12/18/30	CH 18	30 ml	01150641
01/01/12/20/30	CH 20	30 ml	01150670
01/01/12/22/30	CH 22	30 ml	01150687
01/01/12/24/30	CH 24	30 ml	01150693